

resDetect™ Human Interleukin-3 (IL-3) ELISA Kit (Residue Testing)

Catalog Number: RES-A010

Pack Size: 96 tests

IMPORTANT: Please carefully read this manual before performing your experiment.

For Research Use Only. Not For Use in Diagnostic or Therapeutic Procedure

RES10-EN.01

ACTO*

INTENDED USE

resDetectTM Human Interleukin-3 (IL-3) ELISA Kit (Residue Testing) was developed for the detection

and quantitative determination of GMP human IL-3 in cell culture supernates.

It is intended for research use only (RUO).

BACKGROUND

Interleukin 3 is also known as IL-3, is a protein that in humans is encoded by the IL3 gene. Interleukin-3

(IL-3) is an interleukin, a type of biological signal (cytokine) that can improve the body's natural response

to disease as part of the immune system. It acts by binding to the interleukin-3 receptor.

To support the development of CAR-T drugs, ACROBiosystems developed Human Interleukin-3 (IL-3)

ELISA Kit (Residue Testing) with rigorous methodological validation, which is used evaluation the

quality of CAR-T products in drug development and CMC quality control stages. Besides, this kit can

also be used for the quantitative determination of GMP human IL-3 Protein (ACROBiosystems,

Cat#GMP-L03H18) concentrations.

PRINCIPLE OF THE ASSAY

This assay kit is used to measure the levels of human IL-3 by employing a standard sandwich-ELISA

format. The micro-plate in the kit has been pre-coated with Anti-IL-3 Antibody. Firstly, add the standard

samples provided in kit and your samples to the plate, incubate and wash the wells. Then add the Biotin-

Anti-IL-3 Antibody to the plate and form Antibody-antigen-biotinylated antibody complex, incubate and

wash the wells. Next add Streptavidin-HRP to the plate, incubate and wash the wells. At last, load the

substrate into the wells and monitor solution color from blue to yellow. The reaction is stopped by the

addition of a stop solution and the intensity of the absorbance can be measured at 450 nm and 630 nm.

The OD Value reflects the amount of human IL-3 bound.

PRECAUTIONS

1. This kit is for research use only and is not for use in diagnostic or therapeutic applications.

2. The kit is suitable for cell supernatants, serum and plasma samples.

1 / 13





- 3. Do not use reagents past their expiration date.
- 4. Do not mix or substitute reagents with those from other kits or other lot number kits.
- 5. If samples generate values higher than the highest standard, dilute the samples with the appropriate calibrator diluent and repeat the assay. If cell supernatant samples need step dilution, except for the final dilution with diluent, other intermediate dilutions can be in cell culture medium.
- 6. Differences in test results can be caused by a variety of factors, including laboratory operator, pipette usage, plate washing technique, reaction time or temperature, and kit storage.
- 7. This kit is designed to remove or reduce some endogenous interference factors in biological samples, and not all possible influencing factors have been removed.

MATERIALS PROVIDED

Table 1. Materials provided

| Catalog | Components | Size | Format | Storage | | |
|------------|--|-----------------------|--------|-----------------------|-----------------------|--|
| Catalog | Components | Components (96 tests) | | Unopened | Opened | |
| RES010-C01 | Pre-coated Anti-IL-3 Antibody Microplate | 1 plate | Solid | 2-8°C | 2-8°C | |
| RES010-C02 | Human IL-3 Standard | 20 μg | Powder | 2-8°C | -70°C | |
| RES010-C03 | Biotin-Anti-IL-3 Antibody | 20 μg | Powder | 2-8°C | -70°C | |
| RES010-C04 | Streptavidin-HRP | 50 μL | Liquid | 2-8°C, avoid light | 2-8°C, avoid light | |
| RES010-C05 | 10×Washing Buffer | 50 mL | Liquid | 2-8°C | 2-8°C | |
| RES010-C06 | 2×Dilution Buffer | 50 mL | Liquid | 2-8°C | 2-8°C | |
| RES010-C07 | Substrate Solution | 12 mL | Liquid | 2-8°C, avoid light | 2-8°C, avoid light | |
| RES010-C08 | Stop Solution | 7 mL | Liquid | 2-8°C | 2-8°C | |

Note: It is recommended that Streptavidin-HRP be centrifuged briefly before use to deposit liquid from the tube wall or cap to the bottom of the tube.

2 / 13

E-mail: order@acrobiosystems.com

US and Canada: Tel: +1 800-810-0816 Web: http://www.acrobiosystems.com

Asia and Pacific: Tel: +86 400-682-2521

RES10-EN.01

SRORAGE

1. Unopened kit should be stored at 2°C-8°C upon receiving.

2. Find the expiration date on the outside packaging and do not use reagents past their expiration date.

3. The opened kit should be stored per components table. The shelf life is 30 days from the date of opening.

REAGENTS/EQUIPMENT NEEDED BUT NOT SUPPLIED

Single or multi-channel micropipettes and pipette tips: need to meet 10 μL, 300 μL, 1000 μL injection

requirements;

37°C Incubator;

Single or dual wavelength microplate reader with 450 nm and 630 nm filter;

Tubes: 1.5mL,10mL;

Timer;

Reagent bottle;

Deionized or distilled water.

REAGENT PREPARATION

Bring all reagents and samples to room temperature (20°C-25°C) before use. If crystals have formed in

buffer solution, place the sample in an 37°C incubator until the crystals have completely dissolved and

bring the solution back to room temperature before use.

According to Table 2, prepare the provided lyophilized product into a storage solution with ultrapure water,

dissolve at room temperature for 15 to 30 minutes, and mix by gently pipetting, avoiding vigorous shaking

or vertexing. The reconstituted storage solution should be stored at -70°C. It is recommended that the

number of freezing and thawing should not exceed 1 time, the size of the aliquot should not be less than

 $5 \mu g$.

Note: Considering inevitable minor quantitation variations between protein batches, it is also reasonable

to generate the standard curve with specific lot of proteins used for current production for even better

accuracy.

3 / 13



| Table 2. | Preparation | method |
|----------|-------------|--------|
|----------|-------------|--------|

| ID | Components | Size (96 T) | Storage solution concentration. | Reconstituted water Vol. |
|------------|---------------------------|-------------|---------------------------------|--------------------------|
| RES010-C02 | Human IL-3 Standard | 20 μg | 200 μg/mL | 100 μL |
| RES010-C03 | Biotin-Anti-IL-3 Antibody | 20 μg | 200 μg/mL | 100 μL |

RECOMMENDED SAMPLE PREPARATION

1. Working Solution Preparation

1.1 Preparation of 1×Washing Buffer:

Dilute 50 mL 10×Washing Buffer with ultrapure water/deionized water to 500 mL.

1.2 Preparation of 1×Dilution Buffer:

Dilute 50 mL 2×Dilution Buffer with 1×Washing Buffer to 100 mL.

1.3 Preparation of Biotin-Anti-IL-3 Antibody working fluid:

Dilute Biotin-Anti-IL-3 Antibody to 0.05 μg/mL with 1×Dilution Buffer. Please prepare it for one-time use only.

1.4 Preparation of Streptavidin-HRP working fluid:

Dilute Streptavidin-HRP at 1:2000 with 1×Dilution Buffer. The prepared working fluid should avoid light. Please prepare it for one-time use only.

- 1.5 Sample preparation
- a. If the sample to be tested is the cell supernatant, dilute test sample at 1:5 with 1×Dilution Buffer. The volume ratio of sample to diluent is 1:4.
- b. If the sample to be tested is serum or plasma, dilute test sample at 1:5 with 1×Dilution Buffer. The volume ratio of sample to diluent is 1:4.

2. Preparation of Standard curve

The concentration of the reconstituted human IL-3 Calibrator (RES010-C02) is $200\,\mu\text{g/mL}_{2}$, prepare (Std.-0) by diluting 10 μ L the reconstituted human IL-3 Calibrator into 990 μ L Sample Dilution Buffer, mix gently well. Then prepare Std.-1' by diluting 10 μ L Std.-0 into 990 μ L Sample Dilution Buffer. At last,


4 / 13

 US and Canada:
 Tel: +1 800-810-0816
 Web: http://www.acrobiosystems.com

 Asia and Pacific:
 Tel: +86 400-682-2521
 E-mail: order@acrobiosystems.com



prepare the highest concentration of standard curve, Std.-1 (200 pg/mL), by diluting 10 μ L Std.- 1' into 990 μ L Sample Dilution Buffer. Prepare 1:1 serial dilution for the standard curve as follows: Pipette 300 μ L of Sample Dilution Buffer into each tube. Make sure to mix well every time. Sample Dilution Buffer serves as blank.



3. Add Samples

Add 100 µL Calibrator and samples to each well. For blank Control wells, please add 100 µL Dilution Buffer.

Note: It is recommended to set double holes for samples and standard curves to be tested.

4. Incubation

Seal the plate with microplate sealing film and incubate at room temperature for 1.0 hour.

5. Washing

Remove the remaining solution by aspiration, add 300 μ L of 1×Washing Buffer to each well, soak for 10 s, remove any remaining 1×Washing Buffer: by aspirating or decanting, invert the plate and blot it against paper towels. Repeat the wash step above for three times.

6. Add Biotin-Anti-IL-3 Antibody

For all wells, add 100 μL Biotin-Anti-IL-3 Antibody (dilute to 0.05 μg/mL) working solution. Please



prepare it for one-time use only.

7. Incubation

Seal the plate with microplate sealing film and incubate at room temperature for 1.0 hour.

8. Washing

Repeat step 5.

9. Add Streptavidin-HRP

For all wells, add 100 µL Streptavidin-HRP (dilute at 1:2000) working solution. Please prepare it for one-time use only, avoid light.

10. Incubation

Seal the plate with microplate sealing film and incubate at room temperature for 30 min.

11. Washing

Repeat step 5.

12. Substrate Reaction

Add 100 µL Substrate Solution to each well. Seal the plate with microplate sealing film and incubate at room temperature for 20 min, avoid light.

13. Termination

Add 50 µL Stop Solution to each well and tap the plate gently to allow thorough mixing.

Note: The color in the wells should change from blue to yellow.

14. Data Recording

Read the absorbance at 450 nm and 630 nm using UV/Vis microplate spectrophotometer within 10 minutes.

Note: To reduce the background noise, subtract the value read at OD_{450nm} with the value read at $OD_{630 nm}$.

CALCULATION OF RESULTS

- 1. Calculate the mean absorbance for each standard, control and sample and subtract average zero standard optical density (OD).
- 2. The standard curve is plotted with the standard concentration as x-axis and the calibrated absorbance value as y-axis. Four parameters logistic are used to draw the standard curve and calculate the sample

6 / 13

 US and Canada:
 Tel: +1 800-810-0816
 Web: http://www.acrobiosystems.com

 Asia and Pacific:
 Tel: +86 400-682-2521
 E-mail: order@acrobiosystems.com





concentration.


3. Normal range of Standard curve: $R^2 \ge 0.9900$.

4. Detection range: 6.25 pg/mL-200 pg/mL. If the OD value of the sample to be tested is higher than 200 pg/mL, the sample shall be diluted with dilution buffer and assay repeated. If the OD value of the sample to be tested is lower than 6.25 pg/mL, the sample should be reported.

Asia and Pacific:




QUICK GUILD





TYPICAL DATA

For each experiment, a standard curve needs to be set for each micro-plate, and the specific OD value may vary depending on different laboratories, testers, or equipments. The following example data is for reference only. The sample concentration was calculated based on the results of the standard curve.



SENSITIVITY

The minimum detectable concentration of IL-3 is 2.590 pg/mL. The minimum detectable concentration was determined by adding twice standard deviations to the mean optical density value of twenty zero standard replicates and calculating the corresponding concentration.

PRECISION

1. Intra-assay Precision

Three samples of known concentration were tested ten times on one plate to assess intra-assay precision.

2. Inter-assay Precision

Three samples of known concentration were tested in three separate assays to assess inter-assay precision.

| | Intra-assay Precision | | | Inter-assay Precision | | |
|--------------|-----------------------|--------|--------|-----------------------|--------|--------|
| Sample | 1 | 2 | 3 | 1 | 2 | 3 |
| n | 10 | 10 | 10 | 3 | 3 | 3 |
| Mean (pg/mL) | 155.413 | 39.269 | 15.203 | 154.111 | 39.556 | 15.137 |

9 / 13

US and Canada:

Tel: +1 800-810-0816

Web: http://www.acrobiosystems.com

Asia and Pacific: Tel: +86 400-682-2521 E-mail: order@acrobiosystems.com



RES10-EN.01

| SD | 4.016 | 0.966 | 0.479 | 3.518 | 0.370 | 0.358 |
|--------|-------|-------|-------|-------|-------|-------|
| CV (%) | 2.6 | 2.5 | 3.1 | 2.3 | 0.9 | 2.4 |

Note: The example data is for reference only.

RECOVERY

Five parts of blank T cell culture supernatant were added with different concentrations of human IL-3, and the T cell culture supernatant without human IL-3 was used as background to calculate the recovery rate. The range of recovery rate is 92.1-97.4%, and the average recovery is 93.8%.

| Sample Type | Average % Recovery | Range |
|----------------------------------|--------------------|-----------|
| T cell culture supernatant (n=5) | 93.8 | 92.1-97.4 |

LINEARITY

To assess the linearity of the assay, samples spiked with high concentrations of IL-3 were serially diluted with calibrator diluent to produce samples with values within the dynamic range of the assay.

| | | Cell culture medium | Cell culture medium | Serum | |
|---------------|----------------------|---------------------|---------------------|-----------|--|
| | | (DMEM) | (1640) | | |
| 1.2 | Average Recovery (%) | 90.6 | 90.3 | 92.4 | |
| 1:2 | Range (%) | 87.0-94.2 | 87.0-94.4 | 90.1-93.5 | |
| 1.4 | Average Recovery (%) | 89.8 | 96.7 | 88.3 | |
| 1:4 Range (%) | | 83.9-94.4 | 93.8-102.2 | 84.0-90.7 | |
| 1:8 | Average Recovery (%) | 90.3 | 85.5 | 91.2 | |
| 1.8 | Range (%) | 89.0-91.2 | 83.7-86.4 | 89.2-92.6 | |
| 1.16 | Average Recovery (%) | 89.8 | 89.9 | 93.0 | |
| 1:16 | Range (%) | 87.5-92.9 | 87.8-93.5 | 88.7-98.0 | |

Note: The example data is for reference only.

10 / 13

US and Canada: Tel: +1 800-810-0816 Web: http://www.acrobiosystems.com

Asia and Pacific: Tel: +86 400-682-2521 E-mail: order@acrobiosystems.com



SPECIFICITY

This assay recognizes natural and recombinant human IL-3. No cross-reactivity was observed when this kit was used to analyze the following recombinant cytokines.

| Human | | | | | |
|-------|-------------|--------------------|------------|--|--|
| IL-2 | IL-12B | SCF | M-CSF | | |
| IL-4 | IL-15 | Flt-3 Ligand | GM-CSF | | |
| IL-5 | IL-17A | Thrombopoietin-TPO | Anti-CD3 | | |
| IL-6 | IL-18 | TGF-beta 1 | Anti-CD28 | | |
| IL-7 | TNF-alpha | VEGF165 | Anti-CD137 | | |
| IL-8 | IFN-gamma | L1R | | | |
| IL-10 | IFN-alpha 1 | BMP-2 | | | |
| IL-11 | FGF basic | G-CSF | | | |

INTERFERING SUBSTANCES

Verify potential matrix effects by adding different levels of DMSO and HSA to the diluted buffer.

| Additive | Tolerated concentration |
|----------|-------------------------|
| DMSO | 10% |
| HSA | 5% |

CALIBRATION

This immunoassay is calibrated against a highly purified *E. coli*-expressed recombinant human IL-3(91/510). Reference Reagent is calibrated by NIBSC/WHO in April 2013.

NIBSC/WHO (91/510) approximate value (U/mL) = $0.003 \times \text{Human IL-3 value (pg/mL)}$.

11 / 13

US and Canada:

Tel: +1 800-810-0816

E-mail: order@acrobiosystems.com



PLATE LAYOUT

| | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 | 10 | 11 | 12 |
|---|-------|-------|--|--|------|--|---------|----|----|--|----|----|
| Α | Std1 | Std1 | () | () | | | | () | () | () | () | |
| В | Std2 | Std2 | | | | | | | () | | | |
| С | Std3 | Std3 | $\left(\begin{array}{c} \cdots \end{array}\right)$ | \bigcirc | | | | | () | $\left(\begin{array}{c} \dots \end{array} \right)$ | | |
| D | Std4 | Std4 | $\left(\begin{array}{c} \cdots \end{array} \right)$ | $\left(\begin{array}{c} \cdots \end{array} \right)$ | | | () | () | () | $\left(\begin{array}{c} \dots \end{array} \right)$ | () | () |
| E | Std5 | Std5 | () | (| | $\left(\begin{array}{c} \cdots \end{array} \right)$ | <u></u> | () | () | $\left(\begin{array}{c} \dots \end{array} \right)$ | () | () |
| F | Std6 | Std6 | () | () | | | <u></u> | () | () | $\left(\begin{array}{c} \dots \end{array} \right)$ | () | () |
| G | Blank | Blank | () | () | | | | () | () | $\left(\begin{array}{c} \dots \\ \end{array} \right)$ | () | () |
| Н | Blank | Blank | () | () | ···) | () | () | () | () | () | () | () |

Note: Blank is a Blank Dilution Buffer hole.

TROUBLESHOOTING GUIDE

| Problem | Cause | Solution |
|--------------------------|----------------------------------|--------------------------------------|
| Poor standard curve | * Inaccurate pipetting | * Check pipettes |
| Longo CV | * Inaccurate pipetting | * Check pipettes |
| Large CV | * Air bubbles in wells | * Remove bubbles in wells |
| High he alremented | * Plate is insufficiently washed | * Review the manual for proper wash. |
| High background | * Contaminated wash buffer | * Make fresh wash buffer |
| Very low readings across | * Incorrect wavelengths | * Check filters/reader |
| the plate | * Insufficient development time | * Increase development time |





| Samples are reading too high, but standard curve looks fine | * Samples contain cytokine levels above assay range | * Dilute samples and run again |
|---|--|---|
| Drift | * Interrupted assay set-up * Reagents not at room temperature | * Assay set-up should be continuous - have all standards and samples prepared appropriately before commencement of theassay * Ensure that all reagents are at room temperature before pipetting into the wells unless otherwise instructed in the antibody inserts |